

BROWN & WILLIAMSON TOBACCO CORPORATION
RESEARCH SERVICES DEPARTMENT

RESTRICTED

EVALUATION OF THE KAMM HOLDER FOR
EFFECTS ON VENTILATION, TAR AND OTHER
CIGARETTE SMOKE DELIVERIES

J.F. Nall
August 31, 1983

Introduction and Summary of Test Results

A. Purpose of the Test

In early 1982, Philip Morris, Inc. undertook a study purportedly designed to measure ventilation of Barclay cigarettes in a group of human test subjects. The tests, as well as later tests performed by Philip Morris' "independent" testing laboratory, the U.S. Testing Company, were conducted with an experimental device especially constructed for testing Barclays, the "Puff Parameter Analyzer" ("PPA").

On the basis of these tests Philip Morris concluded -- and argued to the FTC -- that Barclay's ventilation is reduced by 30% in "normal" human smoking. The FTC has endorsed this view in its current litigation against Brown & Williamson.

Two explanations have been advanced for the alleged reduction in ventilation during smoking of Barclays:

2023100087

first, that smokers crush Barclay's ventilation channels with their lips, and second, that smokers occlude the ventilation channels, blocking them at the very end of the filter tip with their lips.

Our consultant, Dr. Roger Kamm, has determined that the average pressure exerted by human lips during normal smoking is about 34 torr.* Our consultant, Dr. Louis Fine, has determined, from intra-oral fiber optic viewing of Barclay smokers, that occlusion is more likely to occur at relatively shallow insertion depths.

Based in part on the work of Drs. Fine and Kamm in seeking to devise artificial lips for use in cigarette testing, Dr. Kamm has designed an experimental cigarette holder, fitted with a silicone rubber tubing which, when inflated, replicates the pressure exerted by human lips on a cigarette filter during normal smoking. We have used this holder to perform a series of tests to determine Barclay's ventilation rates and tar, nicotine, and carbon monoxide delivery, at a variety of simulated lip pressures and insertion depths. The pressures and insertion depths used in

* Kamm, Roger D., "PRESSURE MEASURING TESTS," Communication to Brown & Williamson Tobacco Corporation. In Dr. Kamm's studies he measured pressure in terms of "torr". Brown & Williamson's test equipment measured pressure in terms of millimeters of mercury ("mm Hg"). These measurements are precisely equivalent, however, because one torr is equal to one mm Hg.

2023100088

these tests reflect our best understanding of actual human lip pressures and insertion depths in smoking.

The Kamm holder was constructed to permit replication, as closely as possible, of the three parameters considered to be of primary importance in affecting ventilation of Barclay cigarettes: lip pressure, depth of filter insertion in the mouth, and "lip drape". Moreover, the tests we performed with the Kamm holder were designed to present a conservative or "worst case" simulation: lip pressures to 85 and 100 torr, well above the 34 torr average known to exist in human smoking, long insertion length (to 10 mm) to maximize opportunities for channel crushing, and short insertion length (to 4 mm) to maximize opportunities for "lip drape". Insertion depth is defined as the length from the filter end of the cigarette that is inserted in the Kamm holder.

The Kamm holders used in this test were fitted with two different moduli of silicone rubber tubing. The modulus of the tubing determines its relative compliance or stretchiness. By changing the modulus of the tubing in the Kamm holder, it is possible to simulate different degrees of "lip drape." The low modulus tubing is the thinnest walled, most compliant material available, and provides for maximum "lip drape," consistent with the "worst case" design of the test.

2023100089

We performed ventilation and smoke delivery tests on Barclay 85 KS cigarettes with the Kamm holder.

A. Ventilation Tests

For purposes of the ventilation study, 12 cigarettes were tested at each of: (a) 15 pressures (0, 15, 25, 30, 35, 40, 45, 50, 55, 60, 65, 70, 75, 85 and 100 mm of Hg., or torr); (b) 4 insertion depths (4, 6, 8 and 10 mm); and (c) two tubings of different modulus (low and medium) for BARCLAY 85 KS cigarettes. The same cigarette was tested at each level of pressure.

Cigarette ventilation data for individual cigarettes are listed in Appendix A. Averages, standard deviations, and sample sizes appear at the righthand side of Table A1 in each row. A single column (read vertically) within an insertion length is the same cigarette. Differences from initial cigarette ventilation at 0 or atmospheric pressure are similarly summarized in Figure A2. Notably, in the worst possible case -- at the shallowest (4 mm) insertion depth and highest (100 torr) pressure with the low modulus tubing -- the average decrement in ventilation was not 30%, but 14.5%.

At conditions giving the largest decrements in ventilation, the ranges of ventilation values from test cigarette to test cigarette are the greatest. Although the

2023100090

cause for this variation is not completely known, it is most likely the result of variations in filter construction and in the axial orientation of the cigarette in the Kamm holder, or of the greater difficulty of sealing the cigarette in the Kamm holder at the relatively shallow insertion depths.

Figures 3-12 grouped the data to examine trends. In these figures, "*" represents the average decrease in cigarette ventilation from its initial reading. The "H" and "L" denote the high and low readings respectively. Figures 11-12 show average changes in ventilation, overlayed for the different insertion length, for each modulus of tubing.* The significant results of these tests are summarized as follows:

1. Cigarette ventilation is most impaired with decreasing insertion length and increasing pressure. (Figures 11 and 12)

2. For the two shallowest insertion lengths used

* Regression analyses were applied to the data in Figures 11-12. The regression model that is shown below was selected for regressing pressure (X) on the change in ventilation (Y).

$$Y = A + BX^2 \quad (1)$$

The coefficients (A and B) are listed in Figure 13. In addition, the correlation coefficients, r^2 , are provided. A perfect correlation between pressure and change in ventilation yields a correlation coefficient of 1.0. Note that all r^2 are close to 1.0. This indicates the model is a good approximation of the relationship between pressure and change in ventilation.

2023100091

in these tests (4 and 6 mm), ventilation decreases more with the low modulus tubing, which simulates maximum "lip drape," than with the medium modulus tubing. (Compare Figures 11 and 12) At the 4 mm depth with maximum (100 torr) pressure with the low modulus tubing, the lowest average rate of ventilation recorded was 71%. However, even at this 4 mm depth, with maximum applied pressure and low modulus tubing, average ventilation rates as high as 86% were recorded. (Figure A1)

3. At the more moderate 8 mm insertion depth, with pressures of 35 torr -- near the average in human smoking -- the average ventilation rate was 85.82%.

4. For both moduli of tubing, test cigarettes inserted to 8 and 10 mm depths yield about the same rates of ventilation across the range of applied pressures. In other words, there is little difference between the 8 and 10 mm insertion depths. (Figures 11, 12)

5. The amount of variation among test cigarettes, in both ventilation and delivery, is most dramatic at the shallowest insertion depths, with the highest applied pressures, especially with the low modulus tubing. (Figures A1, B1)

B. Smoke Delivery Tests

Samples from the same population of BARCLAY KS cigarettes used in the ventilation study were used for smoking

2023100092

experiments. Twelve ports were smoked for each of 18 combinations (2 moduli x 3 insertion lengths x 3 pressures) of test variables.* Each 12 port sample was randomized across two smoking machines and different runs. Smoke deliveries at 0 torr (atmospheric pressure) for this population of BARCLAY were previously determined using standard FTC smoke procedures (see Figure 14).

Summary statistics for smoke deliveries are shown in Figures 15-17. All of the smoking data are listed in Appendix B. Procedures for measuring smoke delivery are given in Appendix E. Data in Appendix B for tar, nicotine, filter nicotine, filter efficiency and water are presented in columns representing six "replicates." The values given for each replicate represent an average result for ten cigarettes. Data for carbon monoxide and puffs are presented in Appendix B in columns representing 12 replicates. The values given for each replicate represent an average result for five cigarettes.

Graphs of smoke deliveries by tubing are provided in Figures 18-31.

* For purposes of the smoking tests, we abandoned the 10 mm insertion depth, which proved to perform much like the 8 mm insertion depth in the ventilation studies. Also, due to time limitations, only 35 torr, 50 torr and 85 torr pressure variables were employed.

2023100093

All smoking data were evaluated using a technique called analysis of variance ("ANOVA"). ANOVAs are summarized in Appendix C.* Generally, ANOVAs showed that average differences greater than 0.1 mg/cig. of tar, 0.02 mg/cig. of nicotine, and 0.1 mg/cig. of CO are statistically significant at $\alpha = 0.05$ (95% confidence level). Other tests of statistical significance for ancillary smoking data are provided in Appendix C.

The significant results of the smoke delivery tests are summarized as follows:

1. In the very worst case, applying simulated lip pressures (up to 85 torr) well above the known average (34 torr), and the shallowest possible insertion depth (4 mm), with maximum simulated lip drape (low modulus tubing), the average tar delivery recorded was 1.5 mg. per cigarette. (Figure B1)

2. At the highest applied pressure (85 torr), and shallowest insertion depth (4 mm) with the low modulus tubing, the average nicotine delivery recorded was .18 mg./cg. (Figure B1)

* ANOVA is a well-recognized statistical method for determining the isolated effect of each specific test variable, and the effects of each possible combination of test variables. See, e.g., B. Ostle, Statistics in Research (Iowa State Univ. Press, 1963).

2023100094

3. At the more moderate 8 mm insertion depth, with applied pressures near the human average (35 torr), with the low modulus tubing, average tar delivery was 0.6 mg/cig., while average nicotine delivery was .07 mg/cig.

4. The data showed that, over the range of applied pressures and insertion depths with both moduli of tubing, tar delivery did not increase to levels even approaching those indicated by Philip Morris' work with the puff parameter analyzer. With this test apparatus, we determined that pressure of 175 torr -- more than five times the average pressure in human smokers -- is required to reduce ventilation on a Barclay 30% (i.e., to achieve a 60% ventilation rate). With our test apparatus, we smoked 2 samples of 16 ports each at 175 torr, using insertion depths of 6 and 8 mm and the medium modulus tubing. At the 6 mm insertion depth, the maximum tar delivery recorded was 5.4 mg. The minimum tar delivery recorded at the same insertion depth and pressure was 3.2 mg., while the average was 4.1 mg. (Figure 32)

5. CO deliveries are consistent with tar and nicotine and follow the same trends. (see Appendix C and Figure 16 for more detail).

Data shown below summarize results from all experiments:

2023100095

Average Tar Delivery (mg/cig.) and Ventilation Rate (%)
at Different Pressures (mm Hg)

Insertion Length (mm)	<u>0 mm</u>		<u>35 mm</u>		<u>50 mm</u>		<u>85 mm</u>	
	<u>Vent</u>	<u>Tar</u>	<u>Vent</u>	<u>Tar</u>	<u>Vent</u>	<u>Tar</u>	<u>Vent</u>	<u>Tar</u>
<u>Low Modulus</u>								
4	85.5	ND*	83.1	0.7	80.4	0.9	73.6	1.5
6	85.3	ND	84.6	0.8	83.5	0.8	78.8	1.5
8	86.2	ND	85.8	0.6	85.4	0.7	83.4	1.0
<u>Medium Modulus</u>								
4	86.8	ND	85.1	0.8	83.8	0.7	78.2	1.3
6	87.1	ND	86.4	0.6	85.8	0.7	83.0	1.1
8	87.7	ND	87.4	0.6	87.1	0.6	86.0	1.0

Medium Modulus, 175 mm Hg

<u>Insertion Length (mm)</u>	<u>Vent</u>	<u>Tar</u>
6	58.9	4.1
8	64.0	2.9

* ND = not determined. At 0 pressure, it was difficult to seal the test cigarettes in the Kamm holder sufficiently to ensure good results in a test for tar delivery. Accordingly, tar values at 0 (atmospheric) pressure were determined by the FTC method, and are set forth in Figure 14.

2023100096

TEST APPARATUS

1. Kamm Assembly. Figure 1 shows the Kamm assembly. The parts are described as follows:

a) Kamm Holder. The Kamm Holder is machined from a cast acrylic rod. Dimensions are given in the figure. The stainless pressure tubing is pressed into the acrylic wall and sealed with epoxy glue. The pressure tubing is .040" ID HPLC stainless tubing manufactured by Waters Assoc.

b) Kamm Collar. The Kamm Collar is machined from a cast acrylic rod. Dimensions are given in the figure.

c) Kamm Tubing. Dimensions of the Kamm tubing are given in the figure. The tubing is made of silicone rubber and is manufactured by: Bivona Surgical Instrument Inc., 5700 West 23rd Avenue, Gary, Indiana 46404. Specifications of the tubing are as follows:

<u>Characteristic</u>	<u>Low Modulus</u>	<u>Medium Modulus</u>
Silicone Material (General Electric)	SE 5218	SE 845
Durometer (Shore A)	25 \pm 5	40 \pm 5
Tensile strength (psi)	1000	1200
Tear strength (psi)	100	200
Ultimate elongation (%)	800	800
Tensile Modulus		
at 50% (psi)	55	
100% (psi)	80	120
200% (psi)	110	

2023100097

2. Pressure/Ventilation Apparatus. The pressure/ventilation apparatus is shown in Figure 2.

3. Cast Acrylic Rod. 0.625" and 0.875" diameter, AIN Plastics, Inc., 249 E. Sanford Blvd., Mt. Vernon, N.Y. 10550.

4. Pressure Tubing. Stainless steel .040" I.D. for HPLC. Waters Associates, Milford, MA 01757.

5. Smoking Machine. Brewer Automatic Pipetting Machine, Model 60470. Modified to take a 35 cc puff volume of 2 second duration.

6. Pressure Transducer. Model DP103-10, Validyne Engineering Corporation, Northridge, CA 91324

7. Laminar Flow Tubes. 1/4" copper tubing. 1118 mm long. Distance between flow measuring point = 203 mm.

8. Analog to Digital Converter. Data Acquisition System A113. Interactive Structure, Inc., Bala Cynwyd, PA.

9. Microcomputer. Apple II Plus.

10. Gelman Holder, Gelman Sciences, Medical Device Division, 674 S. Wagner Avenue, Ann Arbor, MI 48106.

11. Puff Parameter Analyzer. This device was similar to the device used by Philip Morris in its ventilation tests. The glass envelope is approximately 135 mm long with a uniform diameter of about 90 mm. The opening at the end of the envelope that accommodates the ventilation holder is 14 mm in diameter.

2023100098

12. Ventilation Holder. Grip holder fitting for cigarettes manufactured by Filtrona Instruments, England.

13. Pressure Regulator. (To provide air pressure at smoking machine.) Nullmatic 0-100 inches of water.

Moore Products Co., Spring House, PA.

EXPERIMENTAL PROCEDURE

The pressure/ventilation apparatus shown in Figure 2 was used to measure the effects on cigarette ventilation of pressure applied to the Kamm holder. Referring to that figure, the Kamm holder (C) is assembled by inserting the Kamm tubing through the plastic cylinder. An air tight seal is made by folding the tubing back and over the rim of the plastic holder. O-rings are slid over the grooves in the holder to insure an air tight seal. The holder/tubing assembly is then inserted into the Kamm collar (B). A Gelman holder (A) is inserted into the other end of the Kamm collar. Photographs of the Kamm assembly are shown in Appendix D.

The ventilation envelope (D) is a Filtrona cigarette holder with dental dam wrapped over each end and sealed with O-rings. The ventilation opening is attached to a pressure transducer (F). The ventilation envelope is attached to the puff parameter analyzer, a glass envelope (E), the other end of which is connected to a pressure transducer (F). A cigarette is inserted through the ventilation envelope and

2023100099

into the puff parameter analyzer so that the tipping perforations are positioned inside the ventilation envelope. The filter end of the cigarette is then inserted into the Kamm holder to the designated insertion depth. The entire apparatus is now ready to be used for testing. Photographs of the apparatus are shown in Appendix D.

The testing sequence consisted of the following operation:

1. Take a standard puff (35 cc puff volume, 2.0 second duration) on the unlit cigarette.
2. By use of the calibrated pressure transducers, A/D converters, and Apple computer, record puff duration, puff volume through cigarette tobacco section, puff volume through ventilation holes, total puff volume, and percent ventilation. All tests are carried out on an unlit cigarette. Repeat puffing for a total of three puff cycles.
3. Using attached pressure syringe, apply desired amount of pressure as indicated by mercury manometer. Then repeat Steps 1 and 2 on each cigarette that is tested.
4. Repeat steps 1, 2 and 3, applying pressures at fifteen graduated increments ranging from 0 to 100 mg Hg. (torr).

The pressure transducers used to measure flow are assembled and calibrated as follows:

2023100100

Taps are brought out of each of the laminar flow elements, and a Validyne DP103-10 differential pressure transducer is hooked to each. The taps are about 203 mm apart and are the same on both elements. The differential transducers then measure the pressure drop across the 203 mm length of tubing. Since the flow is laminar, this pressure is linearly proportional to the flow through the system.

The transducers are attached to individual carrier/demodulators, giving the ability to adjust both span and zero independently of the vent and tobacco section sides.

The two transducer electrical outputs are attached to two channels of multi-channel A/D converter which are interfaced to an Apple II microcomputer. The computer software is designed to detect puff activity when the vent side signal reaches an operator selectable level for an operator selectable number of data points. When this trigger occurs, both the vent and the tobacco side output are summed separately for the duration of the puff. These operations are performed in machine code for maximum acquisition speed. The raw numbers are then converted into volume readings in a BASIC program.

2023100101

The calibration procedure involved the following steps:

1. Set the transducer output so that it is compatible with the input requirements of the A/D converter.
 2. Hook the two laminar flow elements in series so that the flow through both is the same.
 3. Attach a 50 cc syringe to the combined laminar flow elements so that when the plunger is pressed inward, the resulting positive flow causes positive electrical outputs from the transducers.
 4. For each calibration, take a number of volume/flow rate readings ranging from 3 cc to 50 cc. Relate these to the output voltage from the transducers using a standard least square fit.
 5. Enter the resulting equation into the BASIC program to convert the transducer output into flow rate readings. We have found that the correlation between puff volume and integrated transducer signal is consistently high (typically >0.97).
 6. Confirm the resulting calibration by inserting a cigarette or standard of known ventilation rate and check calibration by using a one port syringe type smoking machine (Brewer automatic pipette) set to a 35 cc/2 sec. puff.
- Cigarettes used in this study were taken from regular production at the Macon Branch. They represent

2023100102

May 2, 1983 production. The cigarettes were conditioned for at least 24 hours prior to testing.* Insertion length marks were placed on the tipping at 4 mm, 6 mm, 8 mm, and 10 mm from the mouth end. Cigarettes from this population were used for both ventilation experiments and for smoking experiments.

Two standard Phipps & Bird smoking machines and standard B&W smoking procedures were used to measure tar at various holder pressures. Standard methods for tar, carbon monoxide, nicotine, and water determination are given in Appendix E of this report. The combined Kamm assembly and Gelman holder was weighed before and after smoking to obtain TPM values. The smoking machine was modified by bringing compressed air through a precision regulator and then to the machine through small Tygon tubing (1/16" I.D.). Each port was fitted with a pressure source that would allow pressures from 0 to 90 torr. Pressure was measured at the source of the pressure just after the regulators and after the last port connection to be sure that there were no leaks during smoking. Photographs of the smoking machine are shown in Appendix D.

* It is standard procedure in the tobacco industry, as is recognized by the FTC, to condition cigarettes before performing tests like these. Cigarettes are conditioned by allowing them to sit at 75° Fahrenheit and 60% relative humidity for at least 24 hours prior to testing.

2023100103


Conclusion

Philip Morris contends that the cigarette ventilation of BARCLAY is reduced to about 60% because of blockage of the grooves by human smokers and that tar delivery is 3-7 mg/cig.


With this apparatus we experimentally determined that a pressure of 175 torr or mm of Hg (3 psi), far exceeding that exerted by human lips, is necessary to obtain a cigarette ventilation of about 60%. We then smoked two samples of 16 ports each at 175 mm of pressure for two insertion lengths (6 mm and 8 mm) and a medium modulus Kamm tubing. Smoking results and cigarette ventilation are shown in Figure 32. At this extreme pressure, tar deliveries increased to 3-4 mg/cig., and nicotine deliveries increased to 0.4-0.5 mg/cig.

0592n


2023100104




KAMM COLLAR
(CAST ACRYLIC)



KAMM TUBE
(SILICONE RUBBER)




KAMM HOLDER
(CAST ACRYLIC)




PRESSURE TUBE
(16 GA. STAINLESS STEEL TUBING)

NOTE: ALL DIMENSIONS IN MM.

 BROWN & WILLIAMSON TOBACCO CORPORATION	
DWN. CLARK APR 8 DATE 8-12-83 SCALE 1:1	TITLE FIGURE 1 KAMM ASSEMBLY

2023100105

Figure 2. Pressure/Ventilation Apparatus



- A. GELMAN HOLDER
- B. COLLAR FOR KAMM HOLDER
- C. KAMM HOLDER
- D. VENTILATION ENVELOPE
- E. GLASS ENVELOPE FOR CIGARETTE
- F. PRESSURE TRANSDUCER

2023100106

Figure 3

PLOT OF EFFECTS ON CIGARETTE VENTILATION AT
DIFFERENT PRESSURES
NOTE: H = MAXIMUM AND L = MINIMUM
HOLDER-12 INSET_LQ-10 MODULUS-LOW




Figure 4

PLOT OF EFFECTS ON CIGARETTE VENTILATION AT
DIFFERENT PRESSURES
NOTE: H = MAXIMUM AND L = MINIMUM
HOLDER=12 INSET_LQ=10 MODULUS=HED





Figure 5

PLOT OF EFFECTS ON CIGARETTE VENTILATION AT
DIFFERENT PRESSURES
NOTE: H = MAXIMUM AND L = MINIMUM
HOLDER=12 INSET_LQ=8 MODULUS=LOW



○ - Indicates Hidden Observation

Figure 6

PLOT OF EFFECTS ON CIGARETTE VENTILATION AT
DIFFERENT PRESSURES
NOTE: H = MAXIMUM AND L = MINIMUM
HOLDER-12 INSET_LQ-8 MODULUS-MED





Figure 7

PLOT OF EFFECTS ON CIGARETTE VENTILATION AT
DIFFERENT PRESSURES

NOTE: H = MAXIMUM AND L = MINIMUM
HOLDER=12 INSET_LQ=6 MODULUS=LOW



() - Indicates Hidden Observation

Figure 8

PLOT OF EFFECTS ON CIGARETTE VENTILATION AT
DIFFERENT PRESSURES
NOTE: H = MAXIMUM AND L = MINIMUM
HOLDER=12 INSET_LQ=8 MODULUS=MED




Figure 9

PLOT OF EFFECTS ON CIGARETTE VENTILATION AT
DIFFERENT PRESSURES

NOTE: H = MAXIMUM AND L = MINIMUM
HOLDER=12 INSET_LQ=4 MODULUS=MED





Figure 10

PLOT OF EFFECTS ON CIGARETTE VENTILATION AT
DIFFERENT PRESSURES
NOTE: H = MAXIMUM AND L = MINIMUM
HOLDER=12 INSET_LQ=4 MODULUS=LOW




○ - Indicates Hidden Observation

Figure 11

PLOT OF AVERAGE EFFECTS ON CIGARETTE VENTILATION AT
DIFFERENT PRESSURES AND INSERTION LENGTHS

NOTE : 1 = 10MM(12), 4 = 4MM(12), 6 = 6MM(12), 8 = 8MM(12)
WHERE () IS THE HOLDER SIZE IN MM
MODULUS-LOW



○ - Indicates Hidden Observation

Figure 12

PLOT OF AVERAGE EFFECTS ON CIGARETTE VENTILATION AT
DIFFERENT PRESSURES AND INSERTION LENGTHS

NOTE : 1 = 10MM(12), 4 = 4MM(12), 6 = 6MM(12), 8 = 8MM(12)
WHERE () IS THE HOLDER SIZE IN MM
MODULUS-MED




Figure 13

Regression Analyses for Average Change in
Cigarette Ventilation (%) versus Pressure (mm of Hg)

MODEL: $Y = A + BX^2$ where Y is the average change in cigarette ventilation and X is pressure.

<u>Modulus</u>	<u>Insertion Length (mm)</u>	<u>Coefficients</u>		<u>Correlation Coefficient</u>
		<u>A</u>	<u>B(X10⁻³)</u>	<u>r²</u>
Low	4	0.737797	1.535284	0.97
	6	-0.376379	0.9767059	0.99
	8	-0.162277	0.4009607	0.99
	10	-0.252144	0.2407161	0.92
Medium	4	0.038224	1.017139	0.99
	6	0.012829	0.5707385	0.99
	8	-0.054039	0.2507447	0.97
	10	-0.048470	0.2827123	0.96

2023100117

Figure 14
Initial Smoke Deliveries for BARCLAY KS
Nonmenthol Samples

<u>Measurements</u>	<u>\bar{x}</u>	<u>s</u>	<u>n</u>
Tar (mg/cig.)	0.3	0.1	7
Nicotine (mg/cig.)	0.07	0.03	9
CO (mg/cig.)	0.5	0.1	19
Puffs (/cig.)	8.5	0.3	19

2023100118

Figure 15

SUMMARY STATISTICS FOR SMOKE DELIVERIES BY MODULUS,
INSERTION LENGTHS AND PRESSURES WITH 12MM HOLDER

		VARIABLE=TAR(MG/CIG)		
		RESULT		
		MEAN	STD	N
MODULUS				
LOW		0.83	0.38	54.00
MEDIUM		0.83	0.37	54.00
		INSERTION LENGTH		
MEAN	4MM	0.88	0.47	36.00
	6MM	0.81	0.36	36.00
	8MM	0.74	0.22	36.00
		PRESSURE(MM OF HG)		
MEAN	35MM	0.69	0.17	36.00
	50MM	0.72	0.28	36.00
	85MM	1.22	0.38	36.00
MEAN		0.88	0.38	108.00

		VARIABLE=NICOTINE(MG/CIG)		
		RESULT		
		MEAN	STD	N
MODULUS				
LOW		0.12	0.05	54.00
MEDIUM		0.10	0.04	54.00
		INSERTION LENGTH		
MEAN	4MM	0.12	0.05	36.00
	6MM	0.11	0.05	36.00
	8MM	0.08	0.03	36.00
		PRESSURE(MM OF HG)		
MEAN	35MM	0.08	0.02	36.00
	50MM	0.10	0.03	36.00
	85MM	0.15	0.04	36.00
MEAN		0.11	0.04	108.00

2023100119

Figure 16

SUMMARY STATISTICS FOR SMOKE DELIVERIES BY MODULUS,
INSERTION LENGTHS AND PRESSURES WITH 12MM HOLDER

VARIABLE=CO(MG/CIG)

		RESULT		
		MEAN	STD	N
MODULUS				
LOW		0.51	0.38	108.00
MEDIUM		0.48	0.28	107.00
INSERTION LENGTH				
MEAN	4MM	0.56	0.40	72.00
	6MM	0.50	0.32	72.00
	8MM	0.43	0.22	71.00
PRESSURE(MM OF HG)				
MEAN	35MM	0.32	0.14	72.00
	50MM	0.45	0.24	71.00
	85MM	0.71	0.40	72.00
MEAN		0.48	0.32	215.00

VARIABLE=FILTER_NICOTINE(MG/CIG)

		RESULT		
		MEAN	STD	N
MODULUS				
LOW		0.34	0.07	54.00
MEDIUM		0.34	0.06	54.00
INSERTION LENGTH				
MEAN	4MM	0.35	0.08	36.00
	6MM	0.35	0.06	36.00
	8MM	0.32	0.05	36.00
PRESSURE(MM OF HG)				
MEAN	35MM	0.28	0.04	36.00
	60MM	0.34	0.05	36.00
	85MM	0.38	0.06	36.00
MEAN		0.34	0.06	108.00

2023100120

Figure 17

SUMMARY STATISTICS FOR SMOKE DELIVERIES BY MODULUS,
INSERTION LENGTHS AND PRESSURES WITH 12MM HOLDER

		RESULT		
		MEAN	STD	N
MODULUS				
LOW		76.20	4.05	64.00
MEDIUM		77.41	4.42	64.00
MEAN	INSERTION LENGTH			
	4MM	74.88	3.60	36.00
	6MM	76.36	3.25	36.00
	8MM	77.68	3.69	36.00
MEAN	PRESSURE (MM OF HG)			
	35MM	78.47	3.39	36.00
	60MM	77.68	3.90	36.00
	85MM	72.89	3.81	36.00
MEAN		76.31	4.36	108.00

		RESULT		
		MEAN	STD	N
MODULUS				
LOW		0.17	0.06	64.00
MEDIUM		0.18	0.05	64.00
MEAN	INSERTION LENGTH			
	4MM	0.18	0.05	36.00
	6MM	0.17	0.07	36.00
	8MM	0.16	0.05	36.00
MEAN	PRESSURE (MM OF HG)			
	35MM	0.15	0.06	36.00
	60MM	0.19	0.06	36.00
	85MM	0.18	0.06	36.00
MEAN		0.17	0.06	108.00


		RESULT		
		MEAN	STD	N
MODULUS				
LOW		8.29	0.34	108.00
MEDIUM		8.24	0.39	108.00
MEAN	INSERTION LENGTH			
	4MM	8.19	0.37	72.00
	6MM	8.30	0.30	72.00
	8MM	8.31	0.42	72.00
MEAN	PRESSURE (MM OF HG)			
	35MM	8.24	0.34	72.00
	60MM	8.28	0.36	72.00
	85MM	8.28	0.41	72.00
MEAN		8.27	0.37	216.00

2023100121

Figure 18

PLOT OF PRESSURE EFFECTS ON SMOKE DELIVERIES AT DIFFERENT
INSERTION DEPT'S BY TYPE OF MODULUS FOR 12MM HOLDER

NOTE: 4 = 4MM 6 = 6MM AND 8 = 8MM
MODULUS-LOW




2023100122

Figure 19

PLOT OF PRESSURE EFFECTS ON SMOKE DELIVERIES AT DIFFERENT
INSERTION DEPT'S BY TYPE OF MODULUS FOR 12MM HOLDER

NOTE: 4 = 4MM, 6 = 6MM AND 8 = 8MM
MODULUS=MEDIUM




2023100123

Figure 20

PLOT OF PRESSURE EFFECTS ON SMOKE DELIVERIES AT DIFFERENT
INSERTION DEPT'S BY TYPE OF MODULUS FOR 12MM HOLDER

NOTE: 4 = 4MM, 6 = 6MM AND 8 = 8MM
MODULUS=LOW




2023100124

Figure 21

PLOT OF PRESSURE EFFECTS ON SMOKE DELIVERIES AT DIFFERENT
INSERTION DEPT'S BY TYPE OF MODULUS FOR 12MM HOLDER

NOTE: 4 = 4MM, 5 = 6MM AND 8 = 8MM
MODULUS-MEDIUM




2023100125

Figure 22

PLOT OF PRESSURE EFFECTS ON SMOKE DELIVERIES AT DIFFERENT
INSERTION DEPT'S BY TYPE OF MODULUS FOR 12MM HOLDER

NOTE: 4 = 4MM, 6 = 6MM AND 8 = 8MM
MODULUS-LOW




2023100126

Figure 23

PLOT OF PRESSURE EFFECTS ON SMOKE DELIVERIES AT DIFFERENT
INSERTION DEPT'S BY TYPE OF MODULUS FOR 12MM HOLDER

NOTE: 4 = 4MM, 6 = 6MM AND 8 = 8MM
MODULUS-MEDIUM




2023100127

Figure 24

PLOT OF PRESSURE EFFECTS ON SMOKE DELIVERIES AT DIFFERENT
INSERTION DEPT'S BY TYPE OF MODULUS FOR 12MM HOLDER

NOTE: 4 = 4MM, 6 = 6MM AND 8 = 8MM
MODULUS = LOW




2023100128

Figure 25

PLOT OF PRESSURE EFFECTS ON SMOKE DELIVERIES AT DIFFERENT
INSERTION DEPT'S BY TYPE OF MODULUS FOR 12MM HOLDER

NOTE: 4 = 4MM, 5 = 5MM AND 6 = 6MM
MODULUS-MEDIUM




2023100129

Figure 26

PLOT OF PRESSURE EFFECTS ON SMOKE DELIVERIES AT DIFFERENT
INSERTION DEPT'S BY TYPE OF MODULUS FOR 12MM HOLDER

NOTE: 4 = 4MM, 6 = 6MM AND 8 = 8MM
MODULUS=LOW




○ - Indicates Hidden Observation

Figure 27

PLOT OF PRESSURE EFFECTS ON SMOKE DELIVERIES AT DIFFERENT
INSERTION DEPT'S BY TYPE OF MODULUS FOR 12MM HOLDER

NOTE: 4 = 4MM, 5 = 5MM AND 6 = 6MM
MODULUS-MEDIUM



2023100131

Figure 28

PLOT OF PRESSURE EFFECTS ON SMOKE DELIVERIES AT DIFFERENT
INSERTION DEPT'S BY TYPE OF MODULUS FOR 12MM HOLDER

NOTE: 4 = 4MM, 6 = 6MM AND 8 = 8MM
MODULUS=LOW





Figure 29

PLOT OF PRESSURE EFFECTS ON SMOKE DELIVERIES AT DIFFERENT
INSERTION DEPT'S BY TYPE OF MODULUS FOR 12MM HOLDER

NOTE: 4 = 4MM, 6 = 6MM AND 8 = 8MM
MODULUS-MEDIUM




2023100133

Figure 30

PLOT OF PRESSURE EFFECTS ON SMOKE DELIVERIES AT DIFFERENT
INSERTION DEPT'S BY TYPE OF MODULUS FOR 12MM HOLDER

NOTE: 4 = 4MM 8 = 8MM AND 8 = 8MM
MODULUS=LOW




2023100134

Figure 31

PLOT OF PRESSURE EFFECTS ON SMOKE DELIVERIES AT DIFFERENT
INSERTION DEPT'S BY TYPE OF MODULUS FOR 12MM HOLDER

NOTE: 4 = 4MM, 6 = 6MM AND 8 = 8MM
MODULUS-MEDIUM



2023100135

Figure 32

Experimental Determination of Pressure that is Necessary to
Obtain Approximately 60% Cigarette Ventilation and Its
Effect on Smoke Deliveries

Smoke Deliveries at 175 mm (Hg) of Pressure

Modulus	Insertion Length (mm)	Tar (mg/cig)	Nicotine (mg/cig)	CO (mg/cig)	Filter Nicotine (mg/CO)	Puff (/cig.)	Filter Efficiency (%)	H ₂ O (mg/cig)	Pressure (mm of Hg)	Cigarette Ventilation (%)		
Medium	6 mm	3.2	0.38	2.7	0.63	8.4	62.4	0.48	150	51.3	57.6	48.2
		4.0	0.48	3.8	0.73	8.0	60.3	0.63	163	61.1	59.5	
		5.0	0.58	4.7	0.82	7.9	58.6	0.74	175	62.0	61.0	61.9
		3.8	0.47	3.5	0.74	7.9	61.2	0.70	213	61.6		
		3.5	0.43	2.9	0.68	7.6	61.3	0.46	225	65.7		
		5.4	0.60	5.1	0.86	7.5	58.9	0.66	\bar{x} = 174	58.9		
		3.9	0.48	3.6	0.73	8.3	60.8	0.64	s = 26.0	5.33		
		3.9	0.46	3.5	0.67	7.9	59.3	0.62	n = 10	10		
		\bar{x} = 4.1	0.49	3.7	0.73	7.9	60.3	0.62				
		s = 0.74	0.07	0.82	0.08	0.31	1.31	0.10				
		n = 8	8	8	8	8	8	8				
		3.1	0.36	2.2	0.57	8.3	61.3	0.48	175	65.2	68.3	71.2
		4.8	0.56	4.3	0.82	7.9	59.4	0.57	175	47.5	73.9	65.8
		2.2	0.31	1.9	0.62	8.3	66.7	0.26	175	56.1		
Medium	8 mm	3.2	0.36	2.3	0.60	8.4	62.5	0.28				
		2.7	0.89	2.5	0.65	8.0	42.2	0.44				
		2.2	0.43	2.1	0.55	8.3	56.1	0.40				
		3.4	0.42	2.8	0.68	8.1	61.8	0.38	\bar{x} = 175	64.0		
		1.9	0.25	1.5	0.55	8.3	68.8	0.28	s = 0	9.20		
		\bar{x} = 2.9	0.45	2.5	0.63	8.2	59.9	0.39	n = 7	7		
		s = 0.93	0.20	0.84	0.09	0.18	8.16	0.11				
		n = 8	8	8	8	8	8	8				

2023100136

APPENDIX A
CIGARETTE VENTILATION DATA

2023100137

Figure A1

Source: <https://www.industrydocuments.ucsf.edu/docs/grkm0000>

2023100139

LISTING OF VENTILATION FOR INDIVIDUAL CIGARETTES BY HOLDER LENGTHS, INSERTION LENGTHS AND PRESSURES FOR LOW MODULUS

Figure A1 (cont.)

[illegible]

LISTING OF VENTILATION FOR INDIVIDUAL CIGARETTES BY HOLDER
LEAKAGE, INSERTION LENGTHS AND PRESSURES
FOR MEDIUM MODULUS

Figure A1 (cont.)

CIGARETTE	PRESSURE (MM HG)												N
	VEN- TIL- TIL- VEN-	VEN- TIL- TIL- VEN-	VEN- TIL- TIL- VEN-	VEN- TIL- TIL- VEN-	VEN- TIL- TIL- VEN-	VEN- TIL- TIL- VEN-	VEN- TIL- TIL- VEN-	VEN- TIL- TIL- VEN-	VEN- TIL- TIL- VEN-	VEN- TIL- TIL- VEN-	VEN- TIL- TIL- VEN-	VEN- TIL- TIL- VEN-	
1	12	11	10	9	8	7	6	5	4	3	2	1	
2	12	11	10	9	8	7	6	5	4	3	2	1	
3	12	11	10	9	8	7	6	5	4	3	2	1	
4	12	11	10	9	8	7	6	5	4	3	2	1	
5	12	11	10	9	8	7	6	5	4	3	2	1	
6	12	11	10	9	8	7	6	5	4	3	2	1	
7	12	11	10	9	8	7	6	5	4	3	2	1	
8	12	11	10	9	8	7	6	5	4	3	2	1	
9	12	11	10	9	8	7	6	5	4	3	2	1	
10	12	11	10	9	8	7	6	5	4	3	2	1	
11	12	11	10	9	8	7	6	5	4	3	2	1	
12	12	11	10	9	8	7	6	5	4	3	2	1	
13	12	11	10	9	8	7	6	5	4	3	2	1	
14	12	11	10	9	8	7	6	5	4	3	2	1	
15	12	11	10	9	8	7	6	5	4	3	2	1	
16	12	11	10	9	8	7	6	5	4	3	2	1	
17	12	11	10	9	8	7	6	5	4	3	2	1	
18	12	11	10	9	8	7	6	5	4	3	2	1	
19	12	11	10	9	8	7	6	5	4	3	2	1	
20	12	11	10	9	8	7	6	5	4	3	2	1	
21	12	11	10	9	8	7	6	5	4	3	2	1	
22	12	11	10	9	8	7	6	5	4	3	2	1	
23	12	11	10	9	8	7	6	5	4	3	2	1	
24	12	11	10	9	8	7	6	5	4	3	2	1	
25	12	11	10	9	8	7	6	5	4	3	2	1	
26	12	11	10	9	8	7	6	5	4	3	2	1	
27	12	11	10	9	8	7	6	5	4	3	2	1	
28	12	11	10	9	8	7	6	5	4	3	2	1	
29	12	11	10	9	8	7	6	5	4	3	2	1	
30	12	11	10	9	8	7	6	5	4	3	2	1	
31	12	11	10	9	8	7	6	5	4	3	2	1	
32	12	11	10	9	8	7	6	5	4	3	2	1	
33	12	11	10	9	8	7	6	5	4	3	2	1	
34	12	11	10	9	8	7	6	5	4	3	2	1	
35	12	11	10	9	8	7	6	5	4	3	2	1	
36	12	11	10	9	8	7	6	5	4	3	2	1	
37	12	11	10	9	8	7	6	5	4	3	2	1	
38	12	11	10	9	8	7	6	5	4	3	2	1	
39	12	11	10	9	8	7	6	5	4	3	2	1	
40	12	11	10	9	8	7	6	5	4	3	2	1	
41	12	11	10	9	8	7	6	5	4	3	2	1	
42	12	11	10	9	8	7	6	5	4	3	2	1	
43	12	11	10	9	8	7	6	5	4	3	2	1	
44	12	11	10	9	8	7	6	5	4	3	2	1	
45	12	11	10	9	8	7	6	5	4	3	2	1	
46	12	11	10	9	8	7	6	5	4	3	2	1	
47	12	11	10	9	8	7	6	5	4	3	2	1	
48	12	11	10	9	8	7	6	5	4	3	2	1	
49	12	11	10	9	8	7	6	5	4	3	2	1	
50	12	11	10	9	8	7	6	5	4	3	2	1	
51	12	11	10	9	8	7	6	5	4	3	2	1	
52	12	11	10	9	8	7	6	5	4	3	2	1	
53	12	11	10	9	8	7	6	5	4	3	2	1	
54	12	11	10	9	8	7	6	5	4	3	2	1	
55	12	11	10	9	8	7	6	5	4	3	2	1	
56	12	11	10	9	8	7	6	5	4	3	2	1	
57	12	11	10	9	8	7	6	5	4	3	2	1	
58	12	11	10	9	8	7	6	5	4	3	2	1	
59	12	11	10	9	8	7	6	5	4	3	2	1	
60	12	11	10	9	8	7	6	5	4	3	2	1	
61	12	11	10	9	8	7	6	5	4	3	2	1	
62	12	11	10	9	8	7	6	5	4	3	2	1	
63	12	11	10	9	8	7	6	5	4	3	2	1	
64	12	11	10	9	8	7	6	5	4	3	2	1	
65	12	11	10	9	8	7	6	5	4	3	2	1	
66	12	11	10	9	8	7	6	5	4	3	2	1	
67	12	11	10	9	8	7	6	5	4	3	2	1	
68	12	11	10	9	8	7	6	5	4	3	2	1	
69	12	11	10	9	8	7	6	5	4	3	2	1	
70	12	11	10	9	8	7	6	5	4	3	2	1	
71	12	11	10	9	8	7	6	5	4	3	2	1	
72	12	11	10	9	8	7	6	5	4	3	2	1	
73	12	11	10	9	8	7	6	5	4	3	2	1	
74	12	11	10	9	8	7	6	5	4	3	2	1	
75	12	11	10	9	8	7	6	5	4	3	2	1	
76	12	11	10	9	8	7	6	5	4	3	2	1	
77	12	11	10	9	8	7	6	5	4	3	2	1	
78	12	11	10	9	8	7	6	5	4	3	2	1	
79	12	11	10	9	8	7	6	5	4	3	2	1	
80	12	11	10	9	8	7	6	5	4	3	2	1	
81	12	11	10	9	8	7	6	5	4	3	2	1	
82	12	11	10	9	8	7	6	5	4	3	2	1	
83	12	11	10	9	8	7	6	5	4	3	2	1	
84	12	11	10	9	8	7	6	5	4	3	2	1	
85	12	11	10	9	8	7	6	5	4	3	2	1	
86	12	11	10	9	8	7	6	5	4	3	2	1	
87	12	11	10	9	8	7	6	5	4	3	2	1	
88	12	11	10	9	8	7	6	5	4	3	2	1	
89	12	11	10	9	8	7	6	5	4	3	2	1	
90	12	11	10	9	8	7	6	5	4	3	2	1	
91	12	11	10	9	8	7	6	5	4	3	2	1	
92	12	11	10	9	8	7	6	5	4	3	2	1	
93	12	11	10	9	8	7	6	5	4	3	2	1	
94	12	11	10	9	8	7	6	5	4	3	2	1	
95	12	11	10	9	8	7	6	5	4	3	2	1	
96	12	11	10	9	8	7	6	5	4	3	2	1	
97	12	11	10	9	8	7	6	5	4	3	2	1	
98	12	11	10	9	8	7	6	5	4	3	2	1	
99	12	11	10	9	8	7	6	5	4	3	2	1	
100	12	11	10	9	8	7	6	5	4	3	2	1	
101	12	11	10	9	8	7	6	5	4	3	2	1	
102	12	11	10	9	8	7	6	5	4	3	2	1	
103	12	11	10	9	8	7	6	5	4	3	2	1	
104	12	11	10	9	8	7	6	5	4	3	2	1	
105	12	11	10	9	8	7	6	5	4	3	2	1	
106	12	11	10	9	8	7	6	5	4	3	2	1	
107	12	11	10	9	8	7	6	5	4	3	2	1	
108	12	11	10	9	8	7	6	5	4	3	2	1	
109	12	11	10	9	8	7	6	5	4	3	2	1	
110	12	11	10	9	8	7	6	5	4	3	2	1	
111	12	11	10	9	8	7	6	5	4	3	2	1	
112	12	11	10	9	8	7	6	5	4	3	2	1	
113	12	11	10	9	8	7	6	5	4	3	2	1	
114	12	11	10	9	8	7	6	5	4	3	2	1	
115	12	11	10	9	8	7	6	5	4	3	2	1	
116	12	11	10	9	8	7	6	5	4	3	2	1	
117	12	11	10	9	8	7	6	5	4	3	2	1	
118	12	11	10	9	8	7	6	5	4	3	2	1	
119	12	11	10	9	8	7	6	5	4	3	2	1	
120	12	11	10	9	8	7	6	5	4	3	2	1	

2023100140

(CONTINUED)

Figure A1 (cont.)

Source: <https://www.industrydocuments.ucsf.edu/docs/grkm0000>

HOLDBER 42144

LISTING OF CHANGE IN VENTILATION FOR INDIVIDUAL CIGARETTES
BY HOLDER LENGTHS, INSERTION LENGTHS AND PRESSURES
FOR LOW MODULUS
NOTE: CHANGE = INITIAL - OTHERS FOR EACH CIGARETTE

Source: <https://www.industrydocuments.ucsf.edu/docs/qrk0000>

Figure A2 (cont.)

LISTING OF CHANGE IN VENTILATION FOR INDIVIDUAL CIGARETTES
BY HOLDER LENGTHS, INSERTION LENGTHS AND PRESSURES
FOR LOW MODULUS
NOTE: CHANGE = INITIAL - OTHERS FOR EACH CIGARETTE

HOLDER 12MM

		CIGARETTE												AVG.	STD.	N
INSERTION LENGTH	PRESSURE (MM OF HG)	1	2	3	4	5	6	7	8	9	10	11	12			
		CHG IN VEN-T.	CHG IN VEN-T.	CHG IN VEN-T.	CHG IN VEN-T.	CHG IN VEN-T.	CHG IN VEN-T.	CHG IN VEN-T.	CHG IN VEN-T.	CHG IN VEN-T.	CHG IN VEN-T.	CHG IN VEN-T.	CHG IN VEN-T.	%	%	%
8MM	0MM	0.00	0.00	0.00	0.00	0.00	0.00	0.00	0.00	0.00	0.00	0.00	0.00	0.00	0.00	12.00
	15MM	0.30	0.10	0.00	0.00	0.10	0.20	0.00	0.20	-0.50	-0.10	0.10	0.20	0.05	0.21	12.00
	25MM	1.00	0.20	-0.10	0.10	0.20	0.10	0.00	0.20	-0.40	0.00	0.80	0.10	0.11	0.51	12.00
	30MM	1.70	0.20	-0.20	0.10	0.10	0.40	0.10	0.50	-0.20	0.10	1.10	0.20	0.34	0.55	12.00
	35MM	1.70	0.30	0.00	0.20	0.10	0.40	0.20	0.70	-0.30	0.30	0.30	0.20	0.34	0.48	12.00
	40MM	1.80	0.50	-0.10	0.00	0.00	0.50	0.20	1.00	-0.30	0.20	0.40	0.40	0.38	0.58	12.00
	45MM	2.00	0.40	0.00	0.10	0.20	0.50	0.50	1.60	0.00	0.20	0.30	0.50	0.53	0.63	12.00
	50MM	2.50	0.50	0.10	0.30	0.40	0.70	0.90	1.80	0.00	0.30	0.30	1.00	0.74	0.75	12.00
	55MM	2.10	0.50	0.10	0.50	0.40	1.20	1.30	2.60	0.30	0.50	0.40	2.00	0.88	0.83	12.00
	60MM	2.50	0.50	0.20	0.70	0.50	1.50	1.40	2.70	0.50	0.70	0.50	2.40	1.18	0.88	12.00
	65MM	3.00	0.50	0.20	0.80	0.70	1.80	1.70	3.00	0.90	1.10	0.70	2.50	1.42	0.88	12.00
	70MM	3.30	0.50	0.30	0.90	1.70	2.10	1.80	3.50	1.50	1.30	1.00	2.80	1.73	1.03	12.00
	75MM	3.80	0.80	0.20	1.10	1.70	2.40	2.00	3.80	2.20	2.40	1.30	3.10	2.07	1.15	12.00
	85MM	4.50	1.00	0.50	1.50	2.50	3.30	2.30	4.60	4.20	3.20	2.20	3.60	2.78	1.35	12.00
	100MM	5.80	1.40	0.70	2.40		4.30	3.60	5.80	7.40	4.80	3.40	4.40	4.00	1.88	11.00
10MM	0MM	0.00	0.00	0.00	0.00	0.00	0.00	0.00	0.00	0.00	0.00	0.00	0.00	0.00	0.00	12.00
	15MM	0.20	0.00	0.10	0.10	0.00	0.00	0.10	-0.40	0.00	0.10	-0.10	0.00	0.01	0.15	12.00
	25MM	0.20	0.00	0.00	0.10	0.10	-0.10	0.00	-0.10	0.10	-0.10	-0.10	0.30	0.03	0.13	12.00
	30MM	0.20	-0.10	-0.10	0.10	0.10	-0.10	0.00	-0.20	0.10	0.10	-0.10	0.20	0.02	0.13	12.00
	35MM	0.20	-0.10	-0.10	0.30	0.20	0.00	0.20	0.00	0.10	0.30	0.00	0.40	0.13	0.17	12.00
	40MM	0.20	0.00	-0.10	0.40	0.00	0.00	0.40	-0.10	0.10	0.50	-0.10	0.40	0.14	0.23	12.00
	45MM	0.20	0.10	0.00	0.50	0.00	-0.20	0.30	-0.10	0.10	0.70	-0.10	0.50	0.17	0.28	12.00
	50MM	0.30	-0.10	0.00	0.70	0.10	-0.10	0.40	-0.20	0.00	0.80	-0.10	0.70	0.21	0.36	12.00
	55MM	0.30	0.00	-0.10	0.60	0.30	0.00	0.50	-0.20	0.00	1.10	-0.20	0.80	0.28	0.42	12.00
	60MM	0.40	0.00	-0.20	1.10	0.50	0.10	0.80	-0.10	0.20	1.40	0.30	1.00	0.45	0.51	12.00
	65MM	0.40	0.00	-0.10	1.20	0.80	0.20	0.80	-0.10	0.10	1.70	0.30	1.10	0.54	0.58	12.00
	70MM	0.40	0.20	-0.10	1.30	1.50	0.20	1.00	-0.10	0.20	2.30	0.10	1.30	0.68	0.77	12.00
	75MM	0.50	0.10	-0.20	1.70	2.20	0.20	1.30	-0.10	0.20	3.10	0.20	1.40	0.88	1.05	12.00
	85MM	1.50	0.30	-0.10	3.30	2.60	0.70	2.00	-0.10	0.20	7.10	0.30	2.10	1.67	2.05	12.00
	100MM	2.40	1.30	-0.10	4.80	3.50	1.50	3.50	0.10	0.10	9.50	0.40	2.80	2.50	2.75	12.00

2023100143

LISTING OF CHANGE IN VENTILATION FOR INDIVIDUAL CIGARETTES
BY HOLDER LENGTHS, INSERTION LENGTHS AND PRESSURES
FOR MEDIUM MODULUS
NOTE: CHANGE = INITIAL - OTHERS FOR EACH CIGARETTE

Source: <https://www.industrydocuments.ucsf.edu/docs/qrk0000>

Figure A2 (cont.)

LISTING OF CHANGE IN VENTILATION FOR INDIVIDUAL CIGARETTES
BY HOLDER LENGTHS, INSERTION LENGTHS AND PRESSURES
FOR MEDIUM MODULUS
NOTE: CHANGE = INITIAL - OTHERS FOR EACH CIGARETTE

HOLDER 12MM

		CIGARETTE												AVG.	STD.	N
INSERTION LENGTH	PRESSURE (MM OF HG)	1	2	3	4	5	6	7	8	9	10	11	12			
		CHG IN VEN-T.	CHG IN VEN-T.	CHG IN VEN-T.	CHG IN VEN-T.	CHG IN VEN-T.	CHG IN VEN-T.	CHG IN VEN-T.	CHG IN VEN-T.	CHG IN VEN-T.	CHG IN VEN-T.	CHG IN VEN-T.	CHG IN VEN-T.			
		%	%	%	%	%	%	%	%	%	%	%	%			
8MM	0MM	0.00	0.00	0.00	0.00	0.00	0.00	0.00	0.00	0.00	0.00	0.00	0.00	0.00	0.00	12.00
	15MM	0.20	0.00	0.20	0.40	0.00	0.10	0.00	0.10	0.00	0.20	0.20	0.20	0.13	0.12	12.00
	25MM	0.10	0.10	0.30	0.40	0.10	0.10	0.50	0.10	0.20	0.10	0.20	0.40	0.21	0.19	12.00
	30MM	0.30	0.30	0.40	0.60	0.10	0.10	0.30	0.00	0.30	0.10	0.20	0.60	0.27	0.18	12.00
	35MM	0.40	0.30	0.60	0.60	0.00	0.00	0.50	0.00	0.30	0.10	0.30	0.70	0.32	0.25	12.00
	40MM	0.30	0.10	1.00	0.40	0.20	0.00	0.40	0.10	0.40	0.20	0.20	0.30	0.30	0.25	12.00
	45MM	0.30	0.00	1.20	0.60	0.30	0.20	0.40	0.20	0.70	0.50	0.20	0.50	0.43	0.31	12.00
	50MM	0.50	0.00	1.50	0.80	0.40	0.20	0.50	0.30	1.50	0.40	0.20	0.60	0.57	0.48	12.00
	55MM	0.60	0.00	1.30	0.40	0.40	0.30	0.60	0.30	1.60	0.60	0.30	0.80	0.60	0.45	12.00
	60MM	0.60	0.10	1.70	0.20	0.80	0.40	0.70	0.30	1.80	0.70	0.40	1.00	0.72	0.54	12.00
	65MM	0.60	0.30	1.90	0.20	0.60	0.60	0.70	0.40	2.20	1.30	0.60	1.20	0.88	0.64	12.00
	70MM	0.70	0.10	1.90	0.20	0.90	0.80	0.70	0.50	2.60	1.90	0.80	1.40	1.05	0.75	12.00
	75MM	0.80	0.20	2.00	0.40	0.90	1.00	1.00	0.60	3.10	2.30	1.00	1.70	1.25	0.86	12.00
	85MM	1.10	1.10	2.10	0.60	1.20	1.90	1.30	0.60	3.50	2.40	2.00	2.10	1.55	0.85	12.00
10MM	100MM	1.40	2.10	2.30	0.90	1.70	3.60	2.90	0.90	8.20	3.20	3.20	2.80	2.77	1.84	12.00
	0MM	0.00	0.00	0.00	0.00	0.00	0.00	0.00	0.00	0.00	0.00	0.00	0.00	0.00	0.00	12.00
	15MM	0.10	0.00	0.00	0.10	0.10	0.30	0.10	0.50	0.10	0.30	0.20	0.00	0.15	0.15	12.00
	25MM	0.20	0.20	0.00	0.40	0.30	0.40	0.30	0.60	0.30	0.60	0.30	0.20	0.28	0.22	12.00
	30MM	0.10	0.10	0.10	0.30	0.40	0.50	0.50	0.80	0.40	0.60	0.40	0.40	0.37	0.24	12.00
	35MM	0.10	0.20	0.10	0.10	0.20	0.50	0.70	0.90	0.60	0.80	0.40	0.50	0.41	0.31	12.00
	40MM	0.10	0.20	0.10	0.10	0.30	0.60	0.80	1.10	0.60	0.90	0.30	0.20	0.42	0.37	12.00
	45MM	0.10	0.20	0.20	0.10	0.30	0.60	1.00	1.10	0.60	1.00	0.40	0.30	0.47	0.39	12.00
	50MM	0.20	0.20	0.30	0.10	0.50	0.80	1.30	1.10	0.80	1.00	0.60	0.30	0.58	0.42	12.00
	55MM	0.20	0.30	0.40	0.00	0.50	0.80	1.90	1.40	0.90	1.10	0.60	0.40	0.71	0.55	12.00
	60MM	0.20	0.30	0.40	0.00	0.60	0.90	2.00	1.30	1.10	1.40	0.60	0.50	0.77	0.58	12.00
	65MM	0.30	0.60	0.80	0.10	0.70	1.00	2.50	1.70	1.10	1.70	0.80	0.70	1.00	0.67	12.00
	70MM	0.20	0.90	0.80	0.30	0.70	1.10	3.30	1.70	1.10	1.90	0.70	0.50	1.11	0.86	12.00
	75MM	0.20	1.00	1.60	0.90	0.80	1.20	3.60	2.00	1.30	2.30	0.90	0.60	1.37	0.91	12.00
	85MM	0.30	1.20	3.30	1.90	1.70	1.60	4.10	2.50	2.00	3.30	1.20	0.60	1.87	1.15	12.00
	100MM	0.40	1.70	4.90	2.70	2.70	1.90	4.70	3.20	2.90	4.90	2.10	0.80	3.15	2.12	12.00

2023100145

APPENDIX B
SMOKING DATA

2023100146

Figure B1

LISTING OF SMOKE ANALYSES FOR INDIVIDUAL PORTS BY MODULUS,
INSERTION LENGTHS AND PRESSURES WITH 12MM HOLDER

VARIABLE-TAR(MG/CIQ)

MODULUS LOW

INSERTION LENGTH	PRESSURE(MM OF HG)	REPLICATES						AVG	STD	N
		1	2	3	4	5	6			
		RES-ULT	RES-ULT	RES-ULT	RES-ULT	RES-ULT	RES-ULT			
4MM	35MM	0.70	0.70	0.80	1.00	0.40	0.70	0.72	0.19	6.00
	50MM	0.80	0.50	1.60	0.60	0.90	0.70	0.85	0.39	6.00
	85MM	1.20	1.20	2.20	1.10	1.20	1.80	1.45	0.45	6.00
5MM	35MM	0.70	0.80	0.80	0.80	0.80	0.80	0.78	0.04	6.00
	50MM	0.60	0.90	0.60	0.60	0.90	1.20	0.80	0.24	6.00
	85MM	2.10	1.30	1.20	1.20	1.50	1.50	1.47	0.34	6.00
8MM	35MM	0.40	0.70	0.70	0.70	0.80	0.60	0.63	0.15	6.00
	50MM	0.70	0.60	0.70	0.70	0.60	0.90	0.70	0.11	6.00
	85MM	1.00	0.80	1.00	1.00	0.90	1.10	0.88	0.08	6.00

MODULUS MEDIUM

INSERTION LENGTH	PRESSURE(MM OF HG)	REPLICATES						AVG	STD	N
		1	2	3	4	5	6			
		RES-ULT	RES-ULT	RES-ULT	RES-ULT	RES-ULT	RES-ULT			
4MM	35MM	0.80	0.90	0.50	1.10	0.60	1.00	0.82	0.23	6.00
	50MM	1.20	1.00	1.00	0.30	0.30	0.60	0.72	0.40	6.00
	85MM	1.60	2.30	1.30	1.00	0.80	1.10	1.35	0.54	6.00
5MM	35MM	0.50	0.50	0.60	0.80	0.70	0.50	0.60	0.13	6.00
	50MM	0.90	0.80	0.60	0.50	1.00	0.40	0.70	0.24	6.00
	85MM	0.90	1.20	1.40	1.20	1.20	0.90	1.13	0.20	6.00
8MM	35MM	0.80	0.70	0.50	0.60	0.50	0.40	0.58	0.15	6.00
	50MM	0.50	0.20	0.70	0.60	0.80	0.70	0.58	0.21	6.00
	85MM	0.80	0.80	1.10	0.80	1.30	0.80	0.85	0.21	6.00

2023100147

Figure B1 (cont.)

LISTING OF SMOKE ANALYSES FOR INDIVIDUAL PORTS BY MODULUS,
INSERTION LENGTHS AND PRESSURES WITH 12MM HOLDER

VARIABLE-NICOTINE(MG/CIQ)

MODULUS LOW

INSERTION LENGTH	PRESSURE(MM OF HG)	REPLICATES						AVG	STD	N
		1 RES- ULT	2 RES- ULT	3 RES- ULT	4 RES- ULT	5 RES- ULT	6 RES- ULT			
4MM	35MM	0.09	0.09	0.08	0.11	0.07	0.10	0.09	0.01	6.00
	50MM	0.10	0.09	0.19	0.08	0.11	0.09	0.11	0.04	6.00
	85MM	0.14	0.15	0.25	0.13	0.15	0.23	0.18	0.05	6.00
6MM	35MM	0.09	0.09	0.09	0.10	0.09	0.09	0.09	0.00	6.00
	50MM	0.10	0.12	0.08	0.08	0.12	0.14	0.11	0.02	6.00
	85MM	0.27	0.15	0.15	0.18	0.19	0.17	0.18	0.05	6.00
8MM	35MM	0.05	0.06	0.07	0.09	0.11	0.08	0.07	0.02	6.00
	50MM	0.10	0.08	0.09	0.11	0.09	0.11	0.10	0.01	6.00
	85MM	0.13	0.12	0.12	0.15	0.11	0.15	0.13	0.02	6.00

MODULUS MEDIUM

INSERTION LENGTH	PRESSURE(MM OF HG)	REPLICATES						AVG	STD	N
		1 RES- ULT	2 RES- ULT	3 RES- ULT	4 RES- ULT	5 RES- ULT	6 RES- ULT			
4MM	35MM	0.10	0.10	0.06	0.14	0.07	0.11	0.10	0.03	6.00
	50MM	0.14	0.14	0.15	0.07	0.13	0.07	0.12	0.04	6.00
	85MM	0.18	0.25	0.16	0.13	0.09	0.11	0.15	0.05	6.00
6MM	35MM	0.07	0.06	0.07	0.08	0.04	0.05	0.06	0.02	6.00
	50MM	0.10	0.10	0.09	0.07	0.12	0.07	0.09	0.02	6.00
	85MM	0.12	0.13	0.17	0.15	0.14	0.10	0.13	0.02	6.00
8MM	35MM	0.09	0.08	0.07	0.08	0.07	0.06	0.07	0.01	6.00
	50MM	0.07	0.04	0.09	0.09	0.11	0.09	0.08	0.02	6.00
	85MM	0.10	0.11	0.13	0.10	0.12	0.11	0.11	0.01	6.00

2023100148

2023100149

Figure B1 (cont.)

LISTING OF SMOKE ANALYSES FOR INDIVIDUAL PORTS BY MODULUS,
INSERTION LENGTHS AND PRESSURES WITH 12MM HOLDER

100-443887-1000

[illegible][illegible]

Figure B1 (cont.)

LISTING OF SMOKE ANALYSES FOR INDIVIDUAL PORTS BY MODULUS,
INSERTION LENGTHS AND PRESSURES WITH 12MM HOLDER

[VARIABLE=FILTER NICOTINE(MG/CIG)]

MODULUS LOW

INSERTION LENGTH	PRESSURE (MM OF HG)	REPLICATES						AVG	STD	N
		1	2	3	4	5	6			
		RES- ULT	RES- ULT	RES- ULT	RES- ULT	RES- ULT	RES- ULT			
4MM	35MM	0.30	0.29	0.28	0.32	0.22	0.31	0.29	0.04	6.00
	50MM	0.30	0.30	0.46	0.26	0.29	0.28	0.32	0.07	6.00
	85MM	0.35	0.41	0.54	0.44	0.41	0.54	0.45	0.08	6.00
6MM	35MM	0.30	0.32	0.30	0.26	0.33	0.34	0.31	0.03	6.00
	50MM	0.40	0.34	0.33	0.32	0.37	0.38	0.36	0.03	6.00
	85MM	0.42	0.37	0.38	0.34	0.46	0.48	0.41	0.05	6.00
8MM	35MM	0.26	0.26	0.28	0.32	0.28	0.26	0.28	0.02	6.00
	50MM	0.35	0.31	0.33	0.36	0.30	0.30	0.32	0.03	6.00
	85MM	0.40	0.34	0.35	0.40	0.39	0.38	0.37	0.03	6.00

MODULUS MEDIUM

INSERTION LENGTH	PRESSURE (MM OF HG)	REPLICATES						AVG	STD	N
		1	2	3	4	5	6			
		RES- ULT	RES- ULT	RES- ULT	RES- ULT	RES- ULT	RES- ULT			
4MM	35MM	0.30	0.33	0.30	0.40	0.29	0.37	0.33	0.04	6.00
	50MM	0.44	0.39	0.36	0.27	0.32	0.31	0.35	0.06	6.00
	85MM	0.43	0.52	0.48	0.36	0.30	0.40	0.41	0.08	6.00
6MM	35MM	0.27	0.25	0.29	0.28	0.32	0.24	0.27	0.03	6.00
	50MM	0.42	0.40	0.37	0.38	0.36	0.28	0.37	0.05	6.00
	85MM	0.33	0.42	0.41	0.35	0.36	0.33	0.37	0.04	6.00
8MM	35MM	0.32	0.34	0.30	0.27	0.26	0.23	0.28	0.04	6.00
	50MM	0.28	0.31	0.38	0.38	0.36	0.32	0.34	0.04	6.00
	85MM	0.38	0.34	0.38	0.27	0.37	0.30	0.34	0.05	6.00

2023100150

Figure B1 (cont.)

LISTING OF SMOKE ANALYSES FOR INDIVIDUAL PORTS BY MODULUS,
INSERTION LENGTHS AND PRESSURES WITH 12MM HOLDER

[VARIABLE:PUFF(CIG)]

MODULUS LOW

INSERTION LENGTH	PRESSURE (MM OF HG)	REPLICATES												AVG	STD	N
		1 RES- ULT	2 RES- ULT	3 RES- ULT	4 RES- ULT	5 RES- ULT	6 RES- ULT	7 RES- ULT	8 RES- ULT	9 RES- ULT	10 RES- ULT	11 RES- ULT	12 RES- ULT			
4MM	35MM	8.20	8.80	8.00	8.40	8.20	8.20	8.10	8.20	8.00	8.70	7.70	8.40	8.24	0.30	12.00
	50MM	8.60	8.80	8.60	8.20	8.40	8.20	8.80	8.60	8.00	7.40	8.70	8.30	8.37	0.40	12.00
	85MM	8.40	8.80	7.80	8.20	7.60	8.20	7.80	7.60	8.00	8.60	8.40	8.40	8.15	0.39	12.00
6MM	35MM	8.20	8.20	8.40	8.00	8.30	8.10	8.20	8.00	8.70	8.30	7.70	8.30	8.20	0.24	12.00
	50MM	7.90	8.40	8.60	8.00	8.00	8.10	8.80	8.00	8.20	8.40	8.40	7.70	8.29	0.38	12.00
	85MM	8.10	8.40	8.20	8.30	8.60	8.30	8.60	8.10	8.20	7.80	8.70	8.40	8.31	0.25	12.00
8MM	35MM	8.40	8.40	8.60	8.60	8.40	8.10	8.10	8.20	8.90	8.40	8.20	8.40	8.38	0.23	12.00
	50MM	8.20	8.80	8.20	8.40	8.40	8.20	8.40	8.40	8.30	8.70	8.40	8.10	8.37	0.21	12.00
	85MM	8.70	8.10	8.80	8.70	8.20	8.10	8.00	8.30	8.50	8.40	8.80	9.00	8.30	0.57	12.00

MODULUS MEDIUM

INSERTION LENGTH	PRESSURE (MM OF HG)	REPLICATES												AVG	STD	N
		1 RES- ULT	2 RES- ULT	3 RES- ULT	4 RES- ULT	5 RES- ULT	6 RES- ULT	7 RES- ULT	8 RES- ULT	9 RES- ULT	10 RES- ULT	11 RES- ULT	12 RES- ULT			
4MM	35MM	8.20	8.10	8.40	7.70	7.80	8.70	8.50	7.90	8.40	8.40	7.80	7.60	8.12	0.36	12.00
	50MM	7.90	8.00	8.10	8.00	7.70	7.50	8.80	8.50	8.50	8.30	8.50	8.20	8.17	0.37	12.00
	85MM	8.20	8.30	7.90	7.80	7.10	8.30	8.70	8.20	8.10	8.20	8.10	8.20	8.09	0.38	12.00
6MM	35MM	8.30	8.20	8.40	8.20	7.80	8.00	8.20	8.60	8.20	8.40	8.80	8.80	8.33	0.28	12.00
	50MM	8.00	7.60	8.40	8.00	8.50	8.10	7.80	8.20	8.40	8.20	8.40	8.60	8.18	0.30	12.00
	85MM	8.50	8.00	8.50	8.80	8.40	8.80	8.80	8.10	8.80	8.70	8.20	8.40	8.60	0.29	12.00
8MM	35MM	8.80	8.30	8.20	8.50	8.50	8.30	8.00	8.60	7.60	8.70	8.00	8.40	8.16	0.52	12.00
	50MM	8.40	7.50	8.00	8.60	8.00	8.50	8.20	8.40	8.00	8.40	7.80	8.40	8.27	0.48	12.00
	85MM	8.80	8.40	8.00	8.10	8.60	8.70	8.00	8.20	8.20	7.50	8.70	8.10	8.36	0.42	12.00

2023100151

Figure B1 (cont.)

LISTING OF SMOKE ANALYSES FOR INDIVIDUAL PORTS BY MODULUS,
INSERTION LENGTHS AND PRESSURES WITH 12MM HOLDER

VARIABLE-FILTER EFFICIENCY(%)

MODULUS LOW

INSERTION LENGTH	PRESSURE(MM OF HG)	REPLICATES						AVG	STD	N
		1	2	3	4	5	6			
		RES-ULT	RES-ULT	RES-ULT	RES-ULT	RES-ULT	RES-ULT			
4MM	35MM	76.9	76.3	77.8	74.4	75.9	75.6	76.15	1.15	6.00
	50MM	75	76.9	70.8	76.5	72.5	76.3	74.67	2.48	6.00
	85MM	71.4	71.9	68.4	77.2	71.9	70.1	71.82	2.86	6.00
6MM	35MM	76.9	78	76.9	72.2	78.6	79.1	76.95	2.49	6.00
	50MM	80	73.9	80.5	80	76.5	73.1	77.17	3.38	6.00
	85MM	60.9	71.2	71.7	68	70.8	73.8	69.40	4.56	6.00
8MM	35MM	83.9	81.3	80	78	71.8	81.3	79.38	4.18	6.00
	50MM	77.8	78.5	78.6	76.6	78.9	73.2	77.10	2.19	6.00
	85MM	76.6	73.9	74.6	73.7	78	70.8	74.20	2.51	6.00

MODULUS MEDIUM

INSERTION LENGTH	PRESSURE(MM OF HG)	REPLICATES						AVG	STD	N
		1	2	3	4	5	6			
		RES-ULT	RES-ULT	RES-ULT	RES-ULT	RES-ULT	RES-ULT			
4MM	35MM	75	76.7	83.3	74.1	80.6	77.1	77.60	3.60	6.00
	50MM	75.9	73.6	70.6	79.4	71.1	81.6	75.37	4.46	6.00
	85MM	70.5	67.5	74.2	73.6	76.9	78.4	73.50	4.03	6.00
6MM	35MM	79.4	80.6	80.6	75.7	88.9	82.8	81.33	4.38	6.00
	50MM	80.8	80	80.4	84.4	75	80	80.10	3.00	6.00
	85MM	73.3	76.4	70.7	70	72	76.7	73.18	2.84	6.00
8MM	35MM	78	81	81.1	77.1	79.8	79.3	79.22	1.60	6.00
	50MM	80.6	88.6	80.9	80.9	78.6	78	80.83	4.18	6.00
	85MM	79.6	78.8	74.6	73	75.6	73.2	75.23	2.41	6.00

2023100152

Figure B1 (cont.)

LISTING OF SMOKE ANALYSES FOR INDIVIDUAL PORTS BY MODULUS,
INSERTION LENGTHS AND PRESSURES WITH 12MM HOLDER

VARIABLE=H2O(MG/CIQ)

MODULUS LOW

INSERTION LENGTH	PRESSURE(MM OF HG)	REPLICATES						AVG	STD	N
		1 RES- ULT	2 RES- ULT	3 RES- ULT	4 RES- ULT	5 RES- ULT	6 RES- ULT			
4MM	35MM	0.24	0.18	0.08	0.16	0.20	0.17	0.17	0.05	6.00
	50MM	0.16	0.23	0.22	0.16	0.16	0.08	0.17	0.05	6.00
	85MM	0.16	0.23	0.24	0.08	0.15	0.25	0.18	0.07	6.00
6MM	35MM	0.16	0.14	0.12	0.10	0.09	0.13	0.12	0.03	6.00
	50MM	0.27	0.19	0.12	0.16	0.28	0.15	0.19	0.07	6.00
	85MM	0.41	0.20	0.24	0.15	0.18	0.14	0.22	0.10	6.00
8MM	35MM	0.14	0.06	0.09	0.14	0.16	0.06	0.11	0.04	6.00
	50MM	0.20	0.16	0.13	0.28	0.17	0.17	0.18	0.05	6.00
	85MM	0.18	0.14	0.18	0.18	0.13	0.14	0.16	0.02	6.00

MODULUS MEDIUM

INSERTION LENGTH	PRESSURE(MM OF HG)	REPLICATES						AVG	STD	N
		1 RES- ULT	2 RES- ULT	3 RES- ULT	4 RES- ULT	5 RES- ULT	6 RES- ULT			
4MM	35MM	0.15	0.15	0.13	0.21	0.18	0.22	0.17	0.04	6.00
	50MM	0.24	0.29	0.28	0.19	0.18	0.15	0.22	0.05	6.00
	85MM	0.16	0.18	0.18	0.22	0.20	0.16	0.18	0.02	6.00
6MM	35MM	0.16	0.12	0.11	0.21	0.14	0.17	0.15	0.04	6.00
	50MM	0.15	0.26	0.22	0.14	0.30	0.19	0.21	0.06	6.00
	85MM	0.18	0.12	0.16	0.16	0.14	0.12	0.14	0.02	6.00
8MM	35MM	0.15	0.15	0.19	0.16	0.30	0.20	0.19	0.05	6.00
	50MM	0.12	0.16	0.21	0.20	0.24	0.13	0.18	0.05	6.00
	85MM	0.18	0.20	0.12	0.19	0.11	0.17	0.16	0.04	6.00

2023100153

APPENDIX C
ANALYSES OF VARIANCE FOR SMOKING DATA

2023100154

Figure C1

ANALYSES OF VARIANCE FOR SMOKE DELIVERIES

VARIABLE=TAR(MG/CIG)

GENERAL LINEAR MODELS PROCEDURE

DEPENDENT VARIABLE: RESULT

SOURCE	DF	SUM OF SQUARES	MEAN SQUARE	F VALUE	PR > F	R-SQUARE	C.V.
MODEL	17	8.54268519	0.50251089	8.79	0.0001	0.551981	30.9542
ERROR	90	8.65833333	0.07388148		ROOT MSE		RESULT MEAN
CORRECTED TOTAL	107	15.20101852			0.27189537		0.87870370

SOURCE	DF	TYPE III SS	F VALUE	PR > F	DF	TYPE IV SS	F VALUE	PR > F
MODULUS	1	0.30083333	4.07	0.0467	1	0.30083333	4.07	0.0467
INSET LG	2	1.14240741	7.72	0.0008	2	1.14240741	7.72	0.0008
MODULUS*INSET_LG	2	0.13722222	0.93	0.3993	2	0.13722222	0.93	0.3993
PRESSURE	2	6.39574074	43.23	0.0001	2	6.39574074	43.23	0.0001
MODULUS*PRESSURE	2	0.05722222	0.38	0.6804	2	0.05722222	0.38	0.6804
INSET LG*PRESSURE	4	0.37537037	1.27	0.2883	4	0.37537037	1.27	0.2883
MODULUS*INSET_LG*PRESSURE	4	0.13388889	0.45	0.7704	4	0.13388889	0.45	0.7704

DUNCAN'S MULTIPLE RANGE TEST FOR VARIABLE: RESULT
NOTE: THIS TEST CONTROLS THE TYPE I COMPARISONWISE ERROR RATE,
NOT THE EXPERIMENTWISE ERROR RATE.

ALPHA=0.05 DF=90 MSE=.0738815

MEANS WITH THE SAME LETTER ARE NOT SIGNIFICANTLY DIFFERENT.

DUNCAN	GROUPING	MEAN	N	MODULUS
	A	0.83148	54	LOW
	B	0.82593	54	MEDIUM

DUNCAN'S MULTIPLE RANGE TEST FOR VARIABLE: RESULT
NOTE: THIS TEST CONTROLS THE TYPE I COMPARISONWISE ERROR RATE,
NOT THE EXPERIMENTWISE ERROR RATE.

ALPHA=0.05 DF=90 MSE=.0738815

MEANS WITH THE SAME LETTER ARE NOT SIGNIFICANTLY DIFFERENT.

DUNCAN	GROUPING	MEAN	N	PRESSURE
	A	1.2222	36	85
	B	0.7250	36	50
	B	0.6889	36	35

DUNCAN'S MULTIPLE RANGE TEST FOR VARIABLE: RESULT
NOTE: THIS TEST CONTROLS THE TYPE I COMPARISONWISE ERROR RATE,
NOT THE EXPERIMENTWISE ERROR RATE.

ALPHA=0.05 DF=90 MSE=.0738815

MEANS WITH THE SAME LETTER ARE NOT SIGNIFICANTLY DIFFERENT.

DUNCAN	GROUPING	MEAN	N	INSET_LG
	A	0.98333	36	4
	A	0.81389	36	6
	B	0.73889	36	8

2023100155

Figure G2

ANALYSES OF VARIANCE FOR SMOKE DELIVERIES

VARIABLE=NICOTINE(MG/CIQ)

GENERAL LINEAR MODELS PROCEDURE

DEPENDENT VARIABLE: RESULT

SOURCE	DF	SUM OF SQUARES	MEAN SQUARE	F VALUE	PR > F	R-SQUARE	C.V.
MODEL	17	0.11846298	0.00696841	8.08	0.0001	0.604414	26.8378
ERROR	90	0.07753333	0.00086148		ROOT MSE		RESULT MEAN
CORRECTED TOTAL	107	0.19599630			0.02935100		0.11018518

SOURCE	DF	TYPE III SS	F VALUE	PR > F	DF	TYPE IV SS	F VALUE	PR > F
MODULUS	1	0.00592593	6.88	0.0102	1	0.00592593	6.88	0.0102
INSET_LQ	2	0.01572407	9.13	0.0002	2	0.01572407	9.13	0.0002
MODULUS*INSET_LQ	2	0.00331296	1.92	0.1522	2	0.00331296	1.92	0.1522
PRESSURE	2	0.08500741	49.34	0.0001	2	0.08500741	49.34	0.0001
MODULUS*PRESSURE	2	0.00311852	1.81	0.1696	2	0.00311852	1.81	0.1696
INSET_LQ*PRESSURE	4	0.00428148	1.24	0.2887	4	0.00428148	1.24	0.2887
MODULUS*INSET_LQ*PRESSURE	4	0.00109259	0.32	0.8659	4	0.00109259	0.32	0.8659

DUNCAN'S MULTIPLE RANGE TEST FOR VARIABLE: RESULT
NOTE: THIS TEST CONTROLS THE TYPE I COMPARISONWISE ERROR RATE,
NOT THE EXPERIMENTWISE ERROR RATE.

ALPHA=0.05 DF=90 MSE=0.6E-04

MEANS WITH THE SAME LETTER ARE NOT SIGNIFICANTLY DIFFERENT.

DUNCAN	GROUPING	MEAN	N	MODULUS
	A	0.11759	54	LOW
	B	0.10278	54	MEDIUM

DUNCAN'S MULTIPLE RANGE TEST FOR VARIABLE: RESULT
NOTE: THIS TEST CONTROLS THE TYPE I COMPARISONWISE ERROR RATE,
NOT THE EXPERIMENTWISE ERROR RATE.

ALPHA=0.05 DF=90 MSE=0.6E-04

MEANS WITH THE SAME LETTER ARE NOT SIGNIFICANTLY DIFFERENT.

DUNCAN	GROUPING	MEAN	N	PRESSURE
	A	0.14833	36	85
	B	0.10056	36	60
	C	0.08167	36	35

DUNCAN'S MULTIPLE RANGE TEST FOR VARIABLE: RESULT
NOTE: THIS TEST CONTROLS THE TYPE I COMPARISONWISE ERROR RATE,
NOT THE EXPERIMENTWISE ERROR RATE.

ALPHA=0.05 DF=90 MSE=0.6E-04

MEANS WITH THE SAME LETTER ARE NOT SIGNIFICANTLY DIFFERENT.

DUNCAN	GROUPING	MEAN	N	INSET_LQ
	A	0.12417	36	4
	A	0.11167	36	8
	B	0.09472	36	8

2023100156

Figure C3

ANALYSES OF VARIANCE FOR SMOKE DELIVERIES

VARIABLE=CO(MQ/CIG)

GENERAL LINEAR MODELS PROCEDURE

DEPENDENT VARIABLE: RESULT

SOURCE	DF	SUM OF SQUARES	MEAN SQUARE	F VALUE	PR > F	R-SQUARE	C.V.
MODEL	17	7.13813568	0.41984916	5.42	0.0001	0.318807	58.2818
ERROR	197	15.25416567	0.07743232				
CORRECTED TOTAL	214	22.39230233			ROOT MSE 0.27826583		RESULT MEAN 0.48441860

SOURCE	DF	TYPE III SS	F VALUE	PR > F	DF	TYPE IV SS	F VALUE	PR > F
MODULUS	1	0.07743200	1.00	0.3185	1	0.07743200	1.00	0.3185
INSET_LQ	2	0.59702569	3.86	0.0228	2	0.59702569	3.86	0.0228
MODULUS*INSET_LQ	2	0.11224792	0.72	0.4857	2	0.11224792	0.72	0.4857
PRESSURE	3	5.69022569	36.74	0.0001	3	5.69022569	36.74	0.0001
MODULUS*PRESSURE	3	0.19989792	1.29	0.2774	3	0.19989792	1.29	0.2774
INSET_LQ*PRESSURE	4	0.35870806	1.16	0.3306	4	0.35870806	1.16	0.3306
MODULUS*INSET_LQ*PRESSURE	4	0.10423061	0.34	0.8531	4	0.10423061	0.34	0.8531

DUNCAN'S MULTIPLE RANGE TEST FOR VARIABLE: RESULT
NOTE: THIS TEST CONTROLS THE TYPE I COMPARISONWISE ERROR RATE,
NOT THE EXPERIMENTWISE ERROR RATE.

ALPHA=0.05 DF=197 MSE=.0774323

WARNING: CELL SIZES ARE NOT EQUAL.
HARMONIC MEAN OF CELL SIZES=107.488

MEANS WITH THE SAME LETTER ARE NOT SIGNIFICANTLY DIFFERENT.

DUNCAN	GROUPING	MEAN	N	MODULUS
	A	0.51288	108	LOW
	A	0.47570	107	MEDIUM

DUNCAN'S MULTIPLE RANGE TEST FOR VARIABLE: RESULT
NOTE: THIS TEST CONTROLS THE TYPE I COMPARISONWISE ERROR RATE,
NOT THE EXPERIMENTWISE ERROR RATE.

ALPHA=0.05 DF=197 MSE=.0774323

WARNING: CELL SIZES ARE NOT EQUAL.
HARMONIC MEAN OF CELL SIZES=71.6638

MEANS WITH THE SAME LETTER ARE NOT SIGNIFICANTLY DIFFERENT.

DUNCAN	GROUPING	MEAN	N	PRESSURE
	A	0.71111	72	85
	B	0.45070	71	60
	C	0.32083	72	35

DUNCAN'S MULTIPLE RANGE TEST FOR VARIABLE: RESULT
NOTE: THIS TEST CONTROLS THE TYPE I COMPARISONWISE ERROR RATE,
NOT THE EXPERIMENTWISE ERROR RATE.

ALPHA=0.05 DF=197 MSE=.0774323

WARNING: CELL SIZES ARE NOT EQUAL.
HARMONIC MEAN OF CELL SIZES=71.6638

MEANS WITH THE SAME LETTER ARE NOT SIGNIFICANTLY DIFFERENT.

DUNCAN	GROUPING	MEAN	N	INSET_LQ
	A	0.56694	72	4
	A	0.48722	72	8
	B	0.42817	71	8

2023100157

Figure C4

ANALYSES OF VARIANCE FOR SMOKE DELIVERIES

VARIABLE=FILTER_NICOTINE(MG/CIG)

GENERAL LINEAR MODELS PROCEDURE

DEPENDENT VARIABLE: RESULT

SOURCE	DF	SUM OF SQUARES	MEAN SQUARE	F VALUE	PR > F	R-SQUARE	C.V.
MODEL	17	0.24086667	0.01416863	6.25	0.0001	0.541314	13.8827
ERROR	90	0.20410000	0.00226776		ROOT MSE		RESULT MEAN
CORRECTED TOTAL	107	0.44496667			0.04762119		0.34277778

SOURCE	DF	TYPE III SS	F VALUE	PR > F	DF	TYPE IV SS	F VALUE	PR > F
MODULUS	1	0.00030000	0.13	0.7169	1	0.00030000	0.13	0.7169
INSET_LQ	2	0.02106667	4.64	0.0120	2	0.02106667	4.64	0.0120
MODULUS*INSET_LQ	2	0.00535556	1.18	0.3117	2	0.00535556	1.18	0.3117
PRESSURE	2	0.17111667	37.73	0.0001	2	0.17111667	37.73	0.0001
MODULUS*PRESSURE	2	0.01567222	3.46	0.0358	2	0.01567222	3.46	0.0358
INSET_LQ*PRESSURE	4	0.02263333	2.60	0.0484	4	0.02263333	2.60	0.0484
MODULUS*INSET_LQ*PRESSURE	4	0.00472222	0.62	0.7208	4	0.00472222	0.62	0.7208

DUNCAN'S MULTIPLE RANGE TEST FOR VARIABLE: RESULT
NOTE: THIS TEST CONTROLS THE TYPE I COMPARISONWISE ERROR RATE,
NOT THE EXPERIMENTWISE ERROR RATE.

ALPHA=0.05 DF=90 MSE=.0022678

MEANS WITH THE SAME LETTER ARE NOT SIGNIFICANTLY DIFFERENT.

DUNCAN	GROUPING	MEAN	N	MODULUS
	A	0.34444	54	LOW
	A	0.34111	54	MEDIUM

DUNCAN'S MULTIPLE RANGE TEST FOR VARIABLE: RESULT
NOTE: THIS TEST CONTROLS THE TYPE I COMPARISONWISE ERROR RATE,
NOT THE EXPERIMENTWISE ERROR RATE.

ALPHA=0.05 DF=90 MSE=.0022678

MEANS WITH THE SAME LETTER ARE NOT SIGNIFICANTLY DIFFERENT.

DUNCAN	GROUPING	MEAN	N	PRESSURE
	A	0.38167	36	85
	B	0.34250	36	60
	C	0.29417	36	35

DUNCAN'S MULTIPLE RANGE TEST FOR VARIABLE: RESULT
NOTE: THIS TEST CONTROLS THE TYPE I COMPARISONWISE ERROR RATE,
NOT THE EXPERIMENTWISE ERROR RATE.

ALPHA=0.05 DF=90 MSE=.0022678

MEANS WITH THE SAME LETTER ARE NOT SIGNIFICANTLY DIFFERENT.

DUNCAN	GROUPING	MEAN	N	INSET_LQ
	A	0.35722	36	4
	A	0.34722	36	6
	B	0.32389	36	8

2023100158

Figure C5

ANALYSES OF VARIANCE FOR SMOKE DELIVERIES

VARIABLE=PUFF(I/CIG)

GENERAL LINEAR MODELS PROCEDURE

DEPENDENT VARIABLE: RESULT

SOURCE	DF	SUM OF SQUARES	MEAN SQUARE	F VALUE	PR > F	R-SQUARE	C.V.
MODEL	17	2.48541667	0.14678922	1.08	0.3763	0.084783	4.4608
ERROR	198	26.83416667	0.13603114				
CORRECTED TOTAL	215	29.42958333					
					ROOT MSE		RESULT MEAN
					0.36882400		0.26805558

SOURCE	DF	TYPE III SS	F VALUE	PR > F	DF	TYPE IV SS	F VALUE	PR > F
MODULUS	1	0.13004630	0.96	0.3284	1	0.13004630	0.96	0.3284
INSET_LQ	2	0.63194444	2.32	0.1007	2	0.63194444	2.32	0.1007
MODULUS*INSET_LQ	2	0.40898148	1.50	0.2249	2	0.40898148	1.50	0.2249
PRESSURE	2	0.07694444	0.28	0.7540	2	0.07694444	0.28	0.7540
MODULUS*PRESSURE	2	0.38453704	1.41	0.2458	2	0.38453704	1.41	0.2458
INSET_LQ*PRESSURE	4	0.61944444	1.14	0.3398	4	0.61944444	1.14	0.3398
MODULUS*INSET_LQ*PRESSURE	4	0.24351852	0.45	0.7741	4	0.24351852	0.45	0.7741

DUNCAN'S MULTIPLE RANGE TEST FOR VARIABLE: RESULT
NOTE: THIS TEST CONTROLS THE TYPE I COMPARISONWISE ERROR RATE,
NOT THE EXPERIMENTWISE ERROR RATE.

ALPHA=0.05 DF=198 MSE=0.136031

MEANS WITH THE SAME LETTER ARE NOT SIGNIFICANTLY DIFFERENT.

DUNCAN	GROUPING	MEAN	N	MODULUS
	A	8.2928	108	LOW
	A	8.2435	108	MEDIUM

DUNCAN'S MULTIPLE RANGE TEST FOR VARIABLE: RESULT
NOTE: THIS TEST CONTROLS THE TYPE I COMPARISONWISE ERROR RATE,
NOT THE EXPERIMENTWISE ERROR RATE.

ALPHA=0.05 DF=198 MSE=0.136031

MEANS WITH THE SAME LETTER ARE NOT SIGNIFICANTLY DIFFERENT.

DUNCAN	GROUPING	MEAN	N	PRESSURE
	A	8.2847	72	85
	A	8.2778	72	50
	A	8.2417	72	35

DUNCAN'S MULTIPLE RANGE TEST FOR VARIABLE: RESULT
NOTE: THIS TEST CONTROLS THE TYPE I COMPARISONWISE ERROR RATE,
NOT THE EXPERIMENTWISE ERROR RATE.

ALPHA=0.05 DF=198 MSE=0.136031

MEANS WITH THE SAME LETTER ARE NOT SIGNIFICANTLY DIFFERENT.

DUNCAN	GROUPING	MEAN	N	INSET_LQ
	A	8.3097	72	8
	A	8.3028	72	6
	A	8.1917	72	4

2023100159

Figure C6

ANALYSES OF VARIANCE FOR SMOKE DELIVERIES
VARIABLE=FILTER_EFFICIENCY(X)
 GENERAL LINEAR MODELS PROCEDURE

DEPENDENT VARIABLE: RESULT

SOURCE	DF	SUM OF SQUARES	MEAN SQUARE	F VALUE	PR > F	R-SQUARE	C.V.
MODEL	17	1066.15333333	62.71490186	5.83	0.0001	0.624131	4.2878
ERROR	90	967.98333334	10.75437037		ROOT MSE		RESULT MEAN
CORRECTED TOTAL	107	2034.13666667			3.27853813		76.30655558

SOURCE	DF	TYPE III SS	F VALUE	PR > F	DF	TYPE IV SS	F VALUE	PR > F
MODULUS	1	131.12037037	12.19	0.0007	1	131.12037037	12.19	0.0007
INSET_LQ	2	140.69555556	6.64	0.0022	2	140.69555556	6.64	0.0022
MODULUS*INSET_LQ	2	30.44740741	1.42	0.2482	2	30.44740741	1.42	0.2482
PRESSURE	2	645.60000000	30.01	0.0001	2	645.60000000	30.01	0.0001
MODULUS*PRESSURE	2	1.29851852	0.06	0.9415	2	1.29851852	0.06	0.9415
INSET_LQ*PRESSURE	4	88.05444444	2.05	0.0945	4	88.05444444	2.05	0.0945
MODULUS*INSET_LQ*PRESSURE	4	29.03703704	0.67	0.6111	4	29.03703704	0.67	0.6111

DUNCAN'S MULTIPLE RANGE TEST FOR VARIABLE: RESULT
 NOTE: THIS TEST CONTROLS THE TYPE I COMPARISONWISE ERROR RATE,
 NOT THE EXPERIMENTWISE ERROR RATE.

ALPHA=0.05 DF=90 MSE=10.7554

MEANS WITH THE SAME LETTER ARE NOT SIGNIFICANTLY DIFFERENT.

DUNCAN	GROUPING	MEAN	N	MODULUS
	A	77.407	54	MEDIUM
	B	75.204	54	LOW

DUNCAN'S MULTIPLE RANGE TEST FOR VARIABLE: RESULT
 NOTE: THIS TEST CONTROLS THE TYPE I COMPARISONWISE ERROR RATE,
 NOT THE EXPERIMENTWISE ERROR RATE.

ALPHA=0.05 DF=90 MSE=10.7554

MEANS WITH THE SAME LETTER ARE NOT SIGNIFICANTLY DIFFERENT.

DUNCAN	GROUPING	MEAN	N	PRESSURE
	A	78.472	36	35
	A	77.556	36	60
	B	72.889	36	85

DUNCAN'S MULTIPLE RANGE TEST FOR VARIABLE: RESULT
 NOTE: THIS TEST CONTROLS THE TYPE I COMPARISONWISE ERROR RATE,
 NOT THE EXPERIMENTWISE ERROR RATE.

ALPHA=0.05 DF=90 MSE=10.7554

MEANS WITH THE SAME LETTER ARE NOT SIGNIFICANTLY DIFFERENT.

DUNCAN	GROUPING	MEAN	N	INSET_LQ
	A	77.678	36	8
	A	76.356	36	6
	B	74.893	36	4

2023100160

Figure C7

ANALYSES OF VARIANCE FOR SMOKE DELIVERIES

VARIABLE=H2O(MQ/CIG)

GENERAL LINEAR MODELS PROCEDURE

DEPENDENT VARIABLE: RESULT

SOURCE	DF	SUM OF SQUARES	MEAN SQUARE	F VALUE	PR > F	II-SQUARE	C.V.
MODEL	17	0.09250463	0.00544145	2.03	0.0171	0.277302	28.7481
ERROR	90	0.24108333	0.00267870		ROOT MSE		RESULT MEAN
CORRECTED TOTAL	107	0.33358786			0.05175618		0.17398148

SOURCE	DF	TYPE III SS	F VALUE	PR > F	DF	TYPE IV SS	F VALUE	PR > F
MODULUS	1	0.00344537	1.29	0.2598	1	0.00344537	1.29	0.2598
INSET_LQ	2	0.00740185	1.38	0.2564	2	0.00740185	1.38	0.2564
MODULUS*INSET_LQ	2	0.00652407	1.22	0.3007	2	0.00652407	1.22	0.3007
PRESSURE	2	0.02696852	6.03	0.0085	2	0.02696852	6.03	0.0085
MODULUS*PRESSURE	2	0.01857963	3.47	0.0354	2	0.01857963	3.47	0.0354
INSET_LQ*PRESSURE	4	0.00755928	0.71	0.5902	4	0.00755928	0.71	0.5902
MODULUS*INSET_LQ*PRESSURE	4	0.02202583	2.08	0.0933	4	0.02202583	2.08	0.0933

DUNCAN'S MULTIPLE RANGE TEST FOR VARIABLE: RESULT
NOTE: THIS TEST CONTROLS THE TYPE I COMPARISONWISE ERROR RATE,
NOT THE EXPERIMENTWISE ERROR RATE.

ALPHA=0.05 DF=90 MSE=.0026787

MEANS WITH THE SAME LETTER ARE NOT SIGNIFICANTLY DIFFERENT.

DUNCAN	GROUPING	MEAN	N	MODULUS
	A	0.17983	54	MEDIUM
	A	0.16833	54	LOW

DUNCAN'S MULTIPLE RANGE TEST FOR VARIABLE: RESULT
NOTE: THIS TEST CONTROLS THE TYPE I COMPARISONWISE ERROR RATE,
NOT THE EXPERIMENTWISE ERROR RATE.

ALPHA=0.05 DF=90 MSE=.0026787

MEANS WITH THE SAME LETTER ARE NOT SIGNIFICANTLY DIFFERENT.

DUNCAN	GROUPING	MEAN	N	PRESSURE
	A	0.19250	36	50
	A	0.17558	36	85
	B	0.16388	36	35

DUNCAN'S MULTIPLE RANGE TEST FOR VARIABLE: RESULT
NOTE: THIS TEST CONTROLS THE TYPE I COMPARISONWISE ERROR RATE,
NOT THE EXPERIMENTWISE ERROR RATE.

ALPHA=0.05 DF=90 MSE=.0026787

MEANS WITH THE SAME LETTER ARE NOT SIGNIFICANTLY DIFFERENT.


DUNCAN	GROUPING	MEAN	N	INSET_LQ
	A	0.18417	36	4
	A	0.17389	36	8
	A	0.16389	36	8

2023100161


APPENDIX D
PHOTO DOCUMENTATION

2023100162


1a. Kamm holder components:
O-rings, Kamm tube, and
Kamm holder




1b. Partially assembled Kamm
holder




1c. Fully assembled Kamm
holder




2023100163



2a. Assembled Kamm holder,
adapter collar, and
Gelman filter




2b. Fully assembled Kamm
holder and Gelman filter




2c. Assembled Kamm holder
and Gelman filter with
cigarette

2023100164


3a. Filtrona ventilation
sleeve with cigarette




3b. Puff profile analyzer
(PPA) with ventilation
sleeve and Kamm holder in
place



c. Complete ventilation
study apparatus



2023100165



4a. Pressure/vacuum manifold
serving two Phipps &
Bird smoke machines



4b. Side view of Phipps &
Bird smoke machine with
Kamm holders in place



4c. Front view of Phipps &
Bird smoke machine with
Kamm holders in place

2023100166